




SUIT reference protocol for OXPHOS analysis by high-resolution respirometry



Doerrier C, Sumbalova Z, Krumschnabel G, Hiller E, Gnaiger E

OROBOROS INSTRUMENTS
high-resolution respirometry
Schöpfstr 18, A-6020 Innsbruck, Austria
Email: instruments@oroboros.at
www.oroboros.at

Contents

1. SUIT RP1: CI(LPE) substrate control	2
2. SUIT RP2: CETF-CI substrate control	5
3. Harmonization between RP1 and RP2.....	8
4. Reference protocols in biological sample.....	8
4.1. RP in HEK 293T cells.....	8
4.2. RP in brain tissue homogenate.....	9
5. Test experiments	10
5.1. SUIT RP1	10
5.2. SUIT RP2	11
6. Technical details	11
7. Author contributions, publication versions, reference	17
8. Appendix.....	Fehler! Textmarke nicht definiert.
Supplement	17

Summary: We developed a substrate-uncoupler-inhibitor titration (SUIT) protocol with the aim to provide a common reference for comparison of respiratory control in mt-preparations obtained from a large variety of species, tissues and cell types. A SUIT reference protocol (SUIT-RP) is required for establishing a database on comparative mitochondrial physiology. The SUIT-RP is applied in the MitoFit proficiency test with HEK 293T cells. It includes a large number of chemicals used in various specific SUIT protocols, subjecting these to quality control in the MitoFit proficiency test.

The SUIT-RP consists of two harmonized SUIT protocols. SUIT-RP1 starts with linear coupling control, $L - P - E$, with the type N substrates pyruvate and malate (N-pathway to Q: CI-linked), thus separating coupling control and the subsequent linear sequence of pathway control in the ETS state (Figure 1). SUIT-RP2 has a focus on OXPHOS capacity of fatty acid oxidation (FAO_P ; F-pathway) compared to OXPHOS capacity with combined NF-type substrates (NF-pathways to Q: $CI\&FAO_P$). RP2 adds a sequence of pathway control steps to measure maximum OXPHOS and ETS capacity with a NFSGp substrate combination to activate pathways converging at the Q-junction through Complex I (CI), electron-transferring flavoprotein complex (CETF), Complex II (CII), and glycerophosphate dehydrogenase complex (CGpDH) ($PGMSOctGp_E$; NFSGp-pathways; Figure 2). Finally, RP1 and RP2 provide information on the activity of the single enzyme step of Complex IV (CIV) downstream of Q. These SUIT-RP are harmonized (Figure 3) such that they can be statistically evaluated as repeat measurements of cross-linked respiratory states, while additional information is obtained when the two protocols are conducted in parallel. Therefore, RP1 and RP2 are complementary with their focus on specific respiratory coupling and pathway control aspects, extending previous strategies for respirometric OXPHOS analysis.

SUIT reference protocols

RP1: 1PM 2D 3c (3NADH) 4U 5G 6S 7Oct 8Rot 9Gp 10Ama 11Tm 12Azd

RP2: 1D 2Oct 3M 4c (4NADH) 5P 6G 7S 8Gp 9U 10Rot 11Ama 12Tm 13Azd

1. SUIT RP1: N(L-P-E) coupling control

PM + mt: NFSGpTm_1PM 2D 3c 4U 5G 6S 7Oct 8Rot 9Gp 10Ama 11Tm 12Azd

	4U	5G	6S	7Oct	8Rot	9Gp	10Ama	11Tm	12Azd
E									
P	2D 3c								
L	1PM								
	N	N	NS	NFS	S	SGp	ROX	CIV	ROX
CI	+	+	+	+	-	-	-	-	-
CETF	-	-	-	+	-	-	-	-	-
CII	-	-	+	+	+	+	-	-	-
CGpDH	-	-	-	-	-	+	-	-	-

Figure 1. RP1: N(L-P-E) coupling control.

SUIT protocol category: NFSGpCIV

SUIT protocol subcategory: N+NS+NFS+S+SGp+CIV

SUIT protocol acronym:

NFSGpTm_1PM 2D 3c 4U 5G 6S 7Oct 8Rot 9Gp 10Ama 11Tm 12Azd

RP1 spotlights

- N or CI-linked linear coupling control: $L - P - E$, thus separating coupling control ($L-P-E$) and pathway control (in the ETS state).

- Oct is added in the ETS state to avoid uncoupler effect with FA.
- N_E (4U and 5G), NS_E (6S), NFS_E (7Oct), S_E (8Rot) and SGp_E (9Gp) are measured. If $NFS_E \approx NS_E$, then this sequence allows calculation of the additivity index of N- and S- linked ETS capacity (related to supercomplex-channeling). This criterium is tested in step 7Oct, evaluating the effect of Oct on NS_E (6S).
- To compare RP1 with RP2: harmonization between protocols in states SGp_E , (RP1:9Gp, RP2:10Rot), and CIV_E (RP1: 11Tm, RP2: 12Tm).
- If Oct is without effect on NS_E (N=PGM; expected in many types of mt), then additional harmonization between protocols is obtained in states PM_p (RP1:2D) = $PMOct_p$ (RP2:5P).
- Harmonization with many previous SUIT protocols up to step 6Rot.

Limitations

- $NFSGp_E$ is not obtained (substrate combination for maximum ETS capacity), in favour of measuring S_E (8Rot). This reference state has to be calculated using the $NFS_E/NFSGp_E$ (9Gp/10Rot; N=PGM) ratio between RP1 and RP2.

RP1mt

mitochondrial preparation: isolated mitochondria (imt), tissue homogenate (thom), and permeabilized fibers (pfi). See Supplement.

Step	State	Comment
PM		CI-linked substrates are added to the medium before mt (mtprep). The state without added substrates is not well defined, slightly higher than ROX due to the presence of some endogenous substrates (shown by a slight decline of respiration and mt-membrane potential upon inhibition by Rot; Krumschnabel et al 2014).
+mt		Incubation up to 20 min to allow stabilization of flux when using high oxygen or during slow exhaustion of endogenous substrates, to obtain N_L .
1PM	PM_L	N-linked LEAK state.
2D	PM_p	OXPHOS coupling efficiency ($P-L$ or $\approx P$ control factor), $j_{\approx P} = \approx P/P = (P-L)/P = 1-L/P$, is measured in the N-linked pathway state (with a possible contribution by partially activating CII-linked respiration; Sumbalova et al 2016a), with defined coupling sites (CI, CIII, CIV) and at high flux compared to FAO (human skeletal muscle: compare with F_p and N_p ; Gnaiger et al 2015).

3c	PM_{cP}	Quality control test early in the SUIT protocol. In pathological states with c release, early addition of c provides information on the $FCF_c = 1-N/N_c$, and separates the FCF_c from other injuries in the subsequent respiratory states. All subsequent states contain c , which is not explicitly written in the following substrate states.
3NADH		NADH is titrated only in case of a high cytochrome c control factor, $FCF_c > 0.1$, to check for the presence of inverted mitochondrial particles and permeability of the inner mt-membrane. If $FCF_c < 0.1$, then the high intactness of the outer mt-membrane provides an indication for intactness of the inner mt-membrane, and thus NADH does not have to be added.
4U	PM_E	CCCP is titrated stepwise to maximum flux, to evaluate limitation of OXPHOS by the phosphorylation system, expressed as the apparent excess $E-P$ capacity factor ($E-P$ coupling control factor), $j_{EXP} = (E-P)/E = 1-P/E$. If $j_{EXP} > 0$, then the ETS coupling efficiency rather than the OXPHOS coupling efficiency is the proper expression of coupling, $j_{\approx E} = \approx E/E = (E-L)/E = 1-L/E$.
5G	PGM_E	$FCR_G = 1-PM/PGM$, reveals an additive effect of convergent electron flux through NADH (N-linked), with a possible contribution by partially activating S-linked respiration.
6S	$PGMS_E$	$FCF_S = 1-N/NS$. It may be important to check if the uncoupler concentration titrated in the PM substrate state is also sufficient for this substrate state.
7Oct	$PGMSOct_E$	$FCF_F = 1-NS/NFS$. This FCF is low or zero in many mt-types. Then also state PM_p is identical to $PMOct_p$ in RP2, and may thus further link the two protocols in the SUIT reference assay for statistical analysis (protocol harmonization). Inhibition is observed at higher FA concentrations.
8Rot	S_E	$FCF_{NF} = 1-S/NFS$. Rot inhibits CI and FAO simultaneously. In some cases it takes very long, until a steady state is reached after inhibition by Rot. Addition of Gp before Rot would not allow a valid estimation of S-linked capacity (compare RP2).
9Gp	SGp_E	Gp-linked capacity is not measured isolated from S in the SUIT reference protocol. $FCF_{Gp} = 1-S/SGp$. This late addition of Gp is a compromise

		for evaluation of the Gp-linked capacity. Malonic acid does not effectively inhibit CII at S_{50} (competitive inhibition at 50 mM succinate). Little is known about the diagnostic value of this Gp-flux control factor. The substrate Gp is expensive.
10Ama	ROX	Inhibition may take very long, particularly in human muscle fibres (Pesta et al 2011; Lemieux et al 2011). This may make ROX correction questionable, particularly if ROX is high in comparison with the initial LEAK state, N_L .
11Tm	CIV_E	Ascorbate (As) is added before TMPD (Tm). $Tm_{0.5}$ is not saturating CIV, and thus represents a compromise, to prevent a too high chemical O_2 background. Apparent CIV activity may thus be lower than ETS capacity determined in the same run. Inhibitor-threshold titrations would be required.
12Azd	ROX	Cyanide is avoided due to the presence of P, but very high Azd concentrations are required. The oxygen dependence of the chemical O_2 is evaluated by a reoxygenation soon after titration of Azd, and can be semi-automatically performed by using the DatLab background calibration function (Slope).

RP1pce (permeabilized cells). See Supplement.

Step	State	Comment
mt	R	In experiments with intact cells (ce), ROUTINE respiration (R) is measured initially, based on endogenous substrates.
PM		CI-linked substrates are added to the medium after mt (ce).
Dig	PM_L	Digitonin permeabilizes the plasma membrane. CI-linked LEAK state of permeabilized cells (pce).

2. SUIT RP2: F-N and NFSGp pathway control

SUIT protocol category: NFSGpCIV

SUIT protocol subcategory: F+NF+NFS+NFSGp+SGp+CIV

SUIT protocol acronym:

NFSGpTm_1D 2Oct 3M 4c 5P 6G 7S 8Gp 9U 10Rot 11Ama 12Tm 13Azd

D+mt: NFSGpTm_1D 2Oct 3M 4c 5P 6G 7S 8Gp 9U 10Rot 11Ama 12Tm 13Azd

						9U	10Rot	11Ama	12Tm	13Azd
E										
P	1D	2Oct 3M 4c	5P	6G	7S	8Gp				
L										
	ROX	F	NF	NF	NFS	NFSGp	SGp	ROX	CIV	ROX
CI	-	-	+	+	+	+	-	-	-	-
CETF	-	+	+	+	+	+	-	-	-	-
CII	-	-	-	-	+	+	+	-	-	-
CGpDH	-	-	-	-	-	+	+	-	-	-

Figure 2. RP2: F-N and NFSGp pathway control.

RP2 spotlights

- Depletion of endogenous substrates with D (1D; State 2).
- F_P (3M) compared to NF_P (5P).
- The full set of pathways converging at Q, NFSGp, is covered (8Gp and 9U in coupling states P and E), and thus the maximum apparent excess E - P capacity factor, $j_{EXP} = 1 - P/E$, can be calculated.
- Harmonization between protocols RP1 and RP2 in states SGp_E , (RP1:9Gp, RP2:10Rot), and CIV_E (RP1: 11Tm, RP2: 12Tm).
- Harmonization with many previous protocols up to S.
- P/E (8Gp/9U) at high ETS capacity compared to RP1.

Limitations

- S_E is not obtained (but it is obtained in RP1).

RP2mt

Mitochondrial preparation: isolated mitochondria (imt), tissue homogenate (thom) and permeabilized fibers (pfi). See Supplement.

Step	State	Comment
+D		In experiments with mt-preparations (mtprep), ADP is added to the medium before the mt.
+mt		D accelerates the depletion of endogenous substrates.
1D	ROX	Substrate depleted ROX state (State 2; Chance, Williams 1955).
2Oct	Oct	Oct alone does not establish an ETS (and OXPHOS) competent substrate state in many mt-types, since M is required to form oxaloacetate and prevent accumulation of acetyl-Co A by the citrate synthase reaction.
	Oct_P	Stimulation of OXPHOS by Oct alone in the presence of D indicates an obscure mechanism of anaplerosis or the presence of N-substrates in the medium.
3M	$OctM_P$	M is titrated stepwise: M.05; M.1; M2. Note that M alone can support OXPHOS if mt-malic enzyme is active, and thus FAO may be overestimated.

4c		See RP1.
4NADH		See RP1.
5P	PMOct _p	M ₂ is required to reduce flux through CII (minimize inhibition by malonate), such that N-linked OXPHOS capacity can be estimated without high scope of compensation by S-linked respiration (Sumbalova et al 2016a). GM _p includes a higher share of S-linked respiration in comparison with PM _p . $FCF_N = 1-F/NF$, important information on training status or cardiac failure (Pesta et al 2011; Lemieux et al, 2011).
6G	PGMOct _p	The state NF _p is obtained.
7S	PGMSOct _p	The state NFS _p is obtained.
8Gp	PGMSOctGp _p	RP2 focuses on maximum P, evaluating additivity of NFSGp in OXPHOS state. $FCF_{Gp} = 1 - NFS/NFSGp$.
9U	PGMSOctGp _E	CCCP is titrated in the NFSGp state with high ETS capacity, to evaluate limitation of OXPHOS by the phosphorylation system. The apparent excess E-P capacity factor (E-P coupling control factor), $j_{EXP} = (E-P)/E = 1-P/E$, is measured in the state of maximum ETS capacity.
10Rot	SGp _E	This state is not a generally valid estimate of S _E (compare RP1 where S _E is obtained). The state SGp _E is identical in RP1 and RP2, and may thus link the two protocols in the SUIT reference assay for statistical analysis (protocol harmonization).
11Ama	ROX	See RP1.
12Tm	CIV _E	See RP1.
13Azd	ROX	See RP1.

RP2pce Permeabilized cells (pce). See Supplement.

Step	State	Comment
mt	R	See RP1pce.
0Dig		See RP1pce.
1D	ROX	See RP1pce.

3. Harmonization between RP1 and RP2




Figure 3. Harmonization of RP1 and RP2.

4. Reference protocols in mitochondrial models

4.1. SUIT-RP in HEK 293T cells




Figure 4. RP1 and RP2 in permeabilized HEK 293T cells. **(A)** RP1 and **(B)** RP2 performed in parallel in MiR06Cr at 37 °C and normoxic conditions. Oxygen concentration ([μM] blue line) and oxygen flow per cells [$\text{pmol}\cdot\text{s}^{-1}\cdot\text{Mill}^{-1}$] (red line).




Figure 5. Harmonization between RP1 and RP2 in permeabilized HEK 293T cells.

4.2. SUIT-RP in brain tissue homogenate




Figure 6. HRR traces obtained in brain tissue homogenate (thom) from mouse after ischemic damage. **(A)** RP1 and **(B)** RP2 carried out in parallel in MiR06Cr at 37 °C and normoxia. Oxygen concentration ([μM] blue line) and oxygen flux per tissue mass [$\text{pmol}\cdot\text{s}^{-1}\cdot\text{mg}^{-1}$] (red line).




Figure 7. Harmonization between RP1 and RP2 in brain thom after ischemic damage.

5. Test experiments

Test experiments are required to finalize the RP2 for specific applications.

- G 10 mM may not be saturating, and higher concentrations should be checked in a test experiment.
- Gp Different sources of Gp are tested (sn-Glycerol 3-phosphate bis(cyclohexylammonium) salt, Sigma G7886 and sn-Glycerol 3-phosphate lithium salt, Sigma, 94124). Gp (sn-Glycerol 3-phosphate lithium salt, 94124 Sigma) is expensive. In preliminary experiments we did not find differences between both Gp types.

5.1. SUIT RP1

- D D is tested to be saturating in NFSGp_p. 7.5 mM (in pfi) may not be saturating in all cases, and higher concentrations of ADP should be checked.
- Depletion of endogenous substrates with D is possible.
- Oct Oct.5 (0.5 mM) might be generally applicable, but in preliminary experiments a higher concentration (1 mM) should be evaluated to check for saturation of flux.
- +S Step titration from S₁₀ to S₅₀ to test if S₁₀ is saturating NS- and S-linked respiratory capacity. If fluxes with S₅₀ > S₁₀ in NS_E, then S₅₀ is added immediately in the OXPHOS state. If S₅₀ < S₁₀, then it is tested if S₅₀ > S₁₀ in S_E, in which case S₅₀ is only added in S_E. If S₁₀ is

saturation in all states, S_{50} may be tested only occasionally, to exclude a shift in the succinate kinetics (in pathologies, ageing, etc).

Dig Optimal digitonin concentration has to be tested in different types of cells for permeabilization of the plasma membrane.
 » <http://wiki.oroboros.at/index.php/Digitonin>

5.2. SUIT RP2

Oct Oct_{0.5} is tested to be saturating in OXPHOS and not inhibiting or uncoupling (titration of high Oct after M.05 or M.1).

M M_{0.1} is tested to be saturating FAO in OXPHOS without activating N-linked respiration beyond F-linked capacity (HEK: mtME). M should be titrated stepwise (M.05; M.1; M2) in the presence of D, to compare the malate kinetics of F_p and N_p .

6. Technical details

Temperature 37 °C.

Data recording interval: 2 s.

Effective chamber volume: 2 ml

Stirrer speed 750 rpm.

DatLab file The default name of the DatLab file contains the date, Power-O2k number and serial experimental number for each day. Example: 2016-01-17 P1-02.DLD

Events Set an 'Event' in DatLab at the time of titration. Use the abbreviated event name, and add information in the comment.

MiR05+CtlCr Catalase (Ctl) is present in all cells, hence addition of Ctl is considered physiological, even if reoxygenations are not required with H₂O₂. MiR06=MiR05+Ctl

Creatine (Cr) is present in many vertebrate cells, and thus should be added generally. With Cr, lower ADP concentrations are saturating for OXPHOS. It may be argued that it should be replaced in invertebrates (*Drosophila*, *C. rabditis*).

MiR06Cr / O2 Mitochondrial respiration medium, 2 ml in the O2k-chamber, plus 100 µl in the capillary of the stopper (more accurately: 88 µl without meniscus). Increase the oxygen concentration to ~450 µM for pfi experiments. Close the chamber.

mt mt-preparation: isolated mitochondria (imt), permeabilized fibers (pfi), tissue homogenate (thom).

pce	Permeabilized cells.
D	If there is time available (20 min), this period may yield a single point for the instrumental high-O ₂ k background. D may be added just before titrating mt (imt, thom) or before opening the chamber for addition of pfi.
pfi / O ₂	During addition of pfi, the O ₂ concentration drops and should be increased immediately to ~450 μM before closing the O ₂ k chamber.
U	'Slope smoothing' should be reduced to 20 (=20 data points used for calculation of the slope), to evaluate very quickly the stimulation of respiration and the need for additional titration steps of CCCP. If only FCCP (more expensive) is available, this can be used and be fully compared with CCCP titrations (a minimally higher CCCP than FCCP concentration may be required for maximum flux).

Cleaning After the experiment clean the O₂k chambers: 3x water, 1x liver homogenate (20 min), 3x water, 3x EtOH 70% (5 min), 1x EtOH 100% (15 min).



O₂k-cleaning SOP

» http://bioblast.at/index.php/MiPNet19.03_O2k-cleaning_and_ISS



Full version with references

» http://wiki.oroboros.at/index.php/MiPNet21.06_SUIT_reference_assay



PM +mt: NFSGpTm_1PM 2D 3c 4U 5G 6S 7Oct 8Rot 9Gp 10Ama 11Tm 12Azd

	4U	5G	6S	7Oct	8Rot	9Gp	10Ama	11Tm	12Azd
E									
P	2D 3c								
L	1PM								
	N	N	NS	NFS	S	SGp	ROX	CIV	ROX
CI	+	+	+	+	-	-	-	-	-
CETF	-	-	-	+	-	-	-	-	-
CII	-	-	+	+	+	+	-	-	-
CGpDH	-	-	-	-	-	+	-	-	-

Event		Mark name	State	Final conc. 2 ml O2k	Stock [mM]	Comment	Tit. [μl]	A	B
O2k and DatLab file: P___ (A / B) 2016- Operator:									
Sample type: Cohort: Sample code:									
Sample.Subsample number: Unit: Concentration:									
Medium:									
MiR									
O2				~200 μM		~450 μM O ₂ for pfi			
P				5 mM	2000		5		
M	0			2 mM	400		10		
mt									
O2	1PM	PM _L				~450 μM O ₂ for pfi			
D	2D	PM _P		2.5 mM	500	7.5 mM (30 μl) for pfi	10		
c	3c	PM _{Pc}		10 μM	4	NADH only if FCF _c >.1	5		
U	4U	PM _E		Δ0.5 μM	1	CCCP	Δ1		
G	5G	PGM _E		10 mM	2000		10		
S	6S	PGMS _E		50 mM	1000		100		
Oct	7Oct	PGMSOct _E		0.5 mM	100		10		
Rot	8Rot	S _E		0.5 μM	1		1		
Gp	9Gp	SGp _E		10 mM	1000		20		
Ama	10Ama	ROX		2.5 μM	5		1		
As				2 mM	800		5		
Tm				0.5 mM	200	~20 min open, C	5		
C	11Tm	Tm _E				~450 μM O ₂ for pfi			
Azd	12Azd	ROX		≥100 mM	4000	~5 min	100		

D+mt: NFSGpTm_1D 2Oct 3M 4c 5P 6G 7S 8Gp 9U 10Rot 11Ama 12Tm 13Azd

E						9U	10Rot	11Ama	12Tm	13Azd
P	1D	2Oct	3M	4c	5P	6G	7S	8Gp		
L										
	ROX	F	NF	NF	NFS	NFSGp	SGp	ROX	CIV	ROX
CI	-	-	+	+	+	+	-	-	-	-
CETF	-	+	+	+	+	+	-	-	-	-
CII	-	-	-	-	+	+	+	-	-	-
CGpDH	-	-	-	-	-	+	+	-	-	-

O2k and DatLab file: P___ (A / B)		2016-	Operator:					
Sample type:		Cohort:	Sample code:					
Sample.Subsample number:		Unit:	Concentration:					
Medium:								
Event	Mark name	State	Final conc. 2 ml O2k	Stock [mM]	Comment	Tit. [μl]	A	B
MiR								
O2			~200 μM		~450 μM for pfi			
D	0		2.5 mM	500	7.5 mM (30 μl) for pfi	10		
mt								
O2	1D	ROX	~200 μM		~450 μM for pfi			
Oct	2Oct	Oct	0.5 mM	100		10		
M.05	3M.05	Oct _p	0.05 mM	50		2		
M.1	3M.1	Oct _p	0.1 mM	50		2		
M2	3M2	Oct _p	2 mM	400		9.5		
c	4c	Oct _{pc}	10 μM	4	NADH only if $FCF_c > .1$	5		
P	5P	PMOct _p	5 mM	2000		5		
G	6G	PGMOct _p	10 mM	2000		10		
S	7S	PGMSOct _p	50 mM	1000		100		
Gp	8Gp	PGMSOctGp _p	10 mM	1000		20		
U	9U	PGMSOctGp _E	Δ0.5 μM	1	CCCP	Δ1		
Rot	10Rot	SGp _E	0.5 μM	1		1		
Ama	11Ama	ROX	2.5 μM	5		1		
As			2 mM	800		5		
Tm			0.5 mM	200	~20 min open, C	5		
C	12Tm	Tm _E			~450 μM O ₂ for pfi, C			
Azd	13Azd	ROX	≥100 mM	4000	~5 min	100		

mt +PM +Dig: NFSGpTm_1PM 2D 3c 4U 5G 6S 7Oct 8Rot 9Gp 10Ama 11Tm 12Azd

	4U	5G	6S	7Oct	8Rot	9Gp	10Ama	11Tm	12Azd
E									
P	2D 3c								
L	1PM								
	N	N	NS	NFS	S	SGp	ROX	CIV	ROX
CI	+	+	+	+	-	-	-	-	-
CETF	-	-	-	+	-	-	-	-	-
CII	-	-	+	+	+	+	-	-	-
CGpDH	-	-	-	-	-	+	-	-	-

O2k and DatLab file: P___ (A / B)		2016-	Operator:					
Sample type:		Cohort:	Sample code:					
Sample.Subsample number:		Unit:	Concentration:					
Medium:								
Event	Mark name	State	Final conc. 2 ml O2k	Stock [mM]	Comment	Tit. [µl]	A	B
MiR								
O2			~200 µM					
mt	R	R						
P			5 mM	2000		5		
M			2 mM	400		10		
Dig	1PM	PM _L		8.1				
D	2D	PM _P	1 / 2.5 mM	500		4 / 10		
c	3c	PM _{Pc}	10 µM	4	NADH only if FCF _c > .1	5		
U	4U	PM _E	Δ0.5 µM	1	CCCP	Δ1		
G	5G	PGM _E	10 mM	2000		10		
S	6S	PGMS _E	50 mM	1000		100		
Oct	7Oct	PGMSOct _E	0.5 mM	100		10		
Rot	8Rot	S _E	0.5 µM	1		1		
Gp	9Gp	SGp _E	10 mM	1000		20		
Ama	10Ama	ROX	2.5 µM	5		1		
As			2 mM	800		5		
Tm			0.5 mM	200	~20 min open, C	5		
C	11Tm	Tm _E						
Azd	12Azd	ROX	≥100 mM	4000	~5 min	100		

mt + Dig: NFSGpTm_1D 2Oct 3M 4c 5P 6G 7S 8Gp 9U 10Rot 11Ama 12Tm 13Azd

E						9U	10Rot	11Ama	12Tm	13Azd
P	1D	2Oct3M4c	5P	6G	7S	8Gp				
L										
	ROX	F	NF	NF	NFS	NFSGp	SGp	ROX	CIV	ROX
CI	-	-	+	+	+	+	-	-	-	-
CETF	-	+	+	+	+	+	-	-	-	-
CII	-	-	-	-	+	+	+	-	-	-
CGpDH	-	-	-	-	-	+	+	-	-	-

O2k and DatLab file: P___ (A / B)		2016-	Operator:					
Sample type:		Cohort:	Sample code:					
Sample.Subsample number:		Unit:	Concentration:					
Medium:								
Event	Mark name	State	Final conc. 2 ml O2k	Stock [mM]	Comment	Tit. [μl]	A	B
MiR								
O2			~200 μM					
mt	R	R						
Dig	0Dig			8.1				
D	1D	ROX	1 / 2.5 mM	500		4 / 10		
Oct	2Oct	Oct	0.5 mM	100		10		
M.05	3M.05	Oct _p	0.05 mM	50		2		
M.1	3M.1	Oct _p	0.1 mM	50		2		
M2	3M2	Oct _p	2 mM	400		9.5		
c	4c	Oct _{pc}	10 μM	4	NADH only if FCF _c > .1	5		
P	5P	PMOct _p	5 mM	2000		5		
G	6G	PGMOct _p	10 mM	2000		10		
S	7S	PGMSOct _p	50 mM	1000		100		
Gp	8Gp	PGMSOctGp _p	10 mM	1000		20		
U	9U	PGMSOctGp _E	Δ0.5 μM	1	CCCP	Δ1		
Rot	10Rot	SGp _E	0.5 μM	1		1		
Ama	11Ama	ROX	2.5 μM	5		1		
As			2 mM	800		5		
Tm			0.5 mM	200	~20 min open, C	5		
C	11Tm	Tm _E						
Azd	13Azd	ROX	≥100 mM	4000	~10 min	100		

Supplement B

B1. Author contributions and acknowledgements

This communication is a pre-publication prepared by CD and EG. CD, ZS, HE, and GK performed test experiments. EG, CD, ZS, and GK contributed to the concept. GK co-wrote the manuscript. EG and CD edited the final version.

We acknowledge the assistance of Giorgia Lamberti for contributing test experiments and editing some figures, and Verena Laner for the MitoFit project management.



Contribution to the project MitoFit, funded by the Tyrolian Government within the programme K-Regio, Standortagentur Tirol, Austria.

<http://www.mitofit.org/index.php/O2k-MitoFit>



B2. General links

Introduction

» http://wiki.orooboros.at/index.php/Gnaiger_2014_MitoPathways

Respiratory substrate-coupling states

» http://www.bioblast.at/index.php/MitoPedia:Respiratory_substrate-coupling_states

Table of titrations

» http://wiki.orooboros.at/index.php/MiPNet09.12_O2k-Titrations

Definition

» http://www.bioblast.at/index.php/Substrate-uncoupler-inhibitor_titration

Context

» http://www.mitofit.org/index.php/SUIT_protocol_library

Abbreviations

» <http://www.bioblast.at/index.php/MitoPedia>

