

Oroboros O2k-Core Manual






Mitochondrial Physiology Network 19.18(C04):1-17 (2016)
 Version C03: 2016-08-24 DatLab 7 ©2014-2016 Oroboros
 C: http://wiki.orooboros.at/index.php/MiPNet19.18C_DatLab-guide

DatLab-guide

Gnaiger E^{1,2}, Plattner C³, Capek O¹, Gradl L⁴

Oroboros Instruments
 Schöpfstrasse 18, 6020 Innsbruck, Austria
 Email: instruments@orooboros.at
www.orooboros.at




Section	1. File	23	Page
	2. Oroboros O2k	25	
	3. TIP2k	25	
	4. Protocol	25	
Experiment	6. Calibration	27	
	7. Flux/Slope	27	
	8. Graph	28	
	9. Layout	33	
	10. Marks.....	36	
	Supplement: Some features of Windows™	39	

This guide through features of DatLab presents an integral component of high-resolution respirometry. For specific applications of DatLab, see:

- » [MiPNet19.18A](#) O2k-start
- » [MiPNet19.18D](#) DatLab O2k-calibration
- » [MiPNet19.18E](#) DatLab O₂ flux analysis: real-time
- » [MiPNet12.10](#) Titration-Injection microPump, TIP2k
- » [MiPNet17.05](#) O2k-Fluo LED2-Module
- » [MiPNet15.03](#) O2k-MultiSensor-ISE
- » [MiPNet15.05](#) NO-manual



1. File

1.1. Open **Ctrl+O** L  Left mouse click **Open** **DLD file** to open a previously saved **DatLab Data file**. An additional DLD file can be opened only in a



separately started DatLab programme.

1.2. Close Ctrl+F4 Close a DLD file. A window **Save changes to file?** pops up offering the options to close the file after saving the changes, or close the file without saving any modifications on the presently open file.

1.3. Save Ctrl+S When disconnected from the O2k, save any changes made under the identical file name overwriting the previous file. Such changes do not affect the raw data of the experiment, but relate to calibrations, experimental protocol, marks, events, and layout.

Temporary backup files are generated by DatLab in the current user's temp directory, indicated by adding tmp.\$\$\$ to the file name. These files are retained only if the PC has failed during data analysis. During data acquisition, the data are written continuously onto the file, hence backup files are not necessary under these conditions.


1.4. Save and disconnect Ctrl+F4 Stop data acquisition and disconnect from the O2k ([MiPNet19.18A](#)).

1.5. Save as When disconnected from the O2k, save the file under a different file name, optionally in a different directory.

1.6. File search Ctrl+F yields a list of all files labelled by the experimental code in a selected directory (see **Experiment \ Edit F3**). Click on the file name to preview the protocol.

1.7. Delete The decision to delete a file containing no useful data can be made most easily when viewing the traces. Only available when disconnected from the O2k.

1.8. Import DatLab templates can be imported for O2k-setups, graph layouts, mark names, TIP2k setups and marks statistics configurations.



DatLab templates can be copied into the programme subdirectory `\DatLab\DLTemplates` from:

http://www.bioblast.at/index.php/DatLab_templates

1.9. Export Data to text file (*.csv) exports plots and events to a text file for further use in Excel and other programs.

Events to text file (*.csv) exports all information in Events to a text file (*.csv). This file may be used as a protocol, including the comments in the Events.


One channel to DatLab 2 analysis exports the O₂ raw signal to a DatLab 2 file (*.DLR); e.g. for O₂ kinetics in DatLab 2.

1.10. Change user Enter the name to update the user code.

1.11. Manage users Rename or delete users.

1.12. Exit Exit DatLab.

2. Oroboros O2k



2.1. O2k-control **F7** Control the O2k-operation mode. » [MiPNet19.18A O2k-start](#)

2.2. O2k-configuration Select or deselect channels that are not actually used, enter sensor numbers and edit channel labels. » [MiPNet19.18A O2k-start](#)

2.3. Stirrer A (B) on/off **F11** (**F12**) Stirrers in chamber A or B are switched on/off.

2.4. Stirrer test **F9** Stirrers are stopped intermittantly (default: 30 s) for a stirrer test. » [MiPNet06.03](#)

2.5. Illumination on/off **F10** The illumination in both chambers is switched on/off.

2.6. Manage setups Setups can be renamed or deleted.

3. TIP2k

» [MiPNet12.10](#) Titration-Injection microPump, TIP2k



4. Protocol

5. Experiment

4.1. Edit **F3** Information on the experimental protocol can be edited at any time during or after the experiment, and all related results are re-calculated instantaneously with the new parameters. Initially, the **Edit experiment**

window displays information from the last file recorded and saved while connected to the O2k.

Reset to system default to reset values to system default.

Cancel to proceed quickly with the experiment, and edit any time later.

Experimental code Up to 10 digits. The File search function **Ctrl+F** lists all files with identical experimental code within a selected directory.

File recorded by (*read only*) shows the user who recorded the file. While connected to the O2k, the **User code** can be changed by **Change user**.

O2k-serial number (*read only*) automatically recorded.



Chamber

(*read only*) as defined in **O2k \ O2k-configuration**.

Protocol Sample

The following entries are entered separately for the left (A) and right (B) O2k-chamber.

Enter the **protocol name**.

Enter information about sample used in each chamber. No sample is added in O₂ calibration experiments.

- **Sample type, Cohort, Sample code, Sample number, Subsample number**

» http://www.bioblast.at/index.php/Edit_experiment_in_DatLab

- **Sample concentration and amount**

Unit Select a unit to express the concentration or amount of sample in the HRR assay.

Million cells - Flow: cell number

mg - Flux: mg of protein, wet weight or dry weight.

Unit - Flux: units of another marker of sample size.

Concentration Enter the sample concentration (e.g. Million cells/ml, mg *W_w*/ml, mg mt-protein/ml). The corresponding amount of sample is calculated on the basis of the O2k-chamber volume.

Amount Alternatively, enter the sample amount (e.g. biopsy *W_w*) if a known amount of sample is added into the chamber. The corresponding sample concentration is calculated on the basis of the chamber volume.

Medium Name of the incubation medium in the O2k-chamber.

Chamber volume The default is 2.00 ml. It is important to define the actually used effective volume of the O2k-chamber for further calculations of oxygen flux.

Data recording interval [s] (*read only*) is selected in the window **O2k-control** **F7**.

Comments For display and printing in the window **Experimental log**.

4.2. Experimental log **Ctrl+F3** The experimental log is generated automatically with information on O2k-settings and calibrations, the **Edit experiment** window and various events. Time-dependent information can be viewed for chambers A, B or both. A filter is selected for viewing minimum, intermittent (default), or all information. **L** **Preview** to view the protocol, and **Save as PDF** file for quality control.



4.3. Add event **F4**

Events An event is a defined point in time, labeled by a name (1 to 10 characters). The event is shown by a vertical line in the graph (line style can be modified under **Graph\Options**) and the label of the event is shown at the top of the graph. A short comment can be entered to describe the event in detail.

Set events **F4** Real-time: Press **F4** to set an event quickly at the current time of the experiment (e.g. to indicate a manual titration into the chamber). The **Edit event** window pops up after setting a new event. Pressing **F4** defines the time point of the event. Full attention can then be paid to the experiment. Edit the event later.

Ctrl+L **L** Insert an event at any chosen moment of the plotted record of the experiment by placing the cursor anywhere in the graph at the selected time point, press **Ctrl** and click the left mouse button **Ctrl+L**.

Edit event **L** Left click on the name of an existing event to open the **Edit event** window to edit or **Delete event**.

Name Enter an event name of 1 to 10 characters. Short names (e.g. O instead of Open) are recommended.

Comment Further information can be entered into the text field.

Chamber A B Both

Select O2k-chamber A, B or both. The Event will be shown on plots for both or one selected chamber.

6. Calibration

Oxygen, potentiometric and amperometric channels

Select chamber and channel to open **Calibration**.

» [MiPNet19.18D O2k-calibration](#)

7. Flux/Slope


Oxygen, potentiometric and amperometric channels

Select chamber and channel to open **Slope configuration**.

» [MiPNet19.18E O2 flux analysis](#)

8. Graph

L The active graph is selected by a left click into the graph. The active graph is highlighted and indicated by the Oroboros logo.




7.1. Graph \ Add A new graph is added at the bottom of the screen. Select plots for display in the new graph, **Ctrl+F6**.

7.2. Graph \ Delete bottom graph The bottom graph is deleted, which reappears with the same layout by **Add**.

7.3. Graph \ Select plots **Ctrl+F6**

7.3.1. Tab: **Oxygen, O2**



Select graph ▼ Pull down to select one of the displayed graphs.

Tabs

Oxygen O2, Amperometric Amp, Potentiometric pX,
System channels: L[Ⓜ] Left click on a channel type to select plots from this channel.

- Y1, Y2** Select the left (Y₁) or right axis (Y₂) for a plot.
- Use default label** selects the default axis label.
- Use custom label** To define a different axis label, select and enter or edit a custom label.

Channel label..Edit (for Amp and pX channels).

- ▼ **Color** Defines the color of the plot.
- ▼ **Width** Defines the line width of the plot.
- ▼ **Line style** Defines the type of line (solid, dash, dot) of the plot.
- ▼ **Unit** Pull down to select a different unit (for Amp and pX).
- Raw signal** displays the non-calibrated raw signal.

- Flux per volume**
- Normalization**
- Flux control ratio**
- Baseline correction**
- Background correction**


[Further details:](http://www.bioblast.at/index.php/Select_plots_in_DatLab)
http://www.bioblast.at/index.php/Select_plots_in_DatLab

Load layout Select layout category (Reference, All users, User: Name), select Layout type: O2k-Core ▼, O2 & Amp ▼', O2 & pX ▼, Other ▼ (depending on the channel selection in O2k-configuration), and Layout name ▼.

Save layout A layout can be saved with any layout name under the category All users or User: Name.

7.3.2. Tab: System channels

- Block temperature** The continuously measured temperature of the copper block, housing the two glass chambers [°C].
- Barometric pressure** The continuously recorded absolute (local) barometric pressure [kPa].
- Ext temp** Signal from external PT1000 temperature sensor (optional), "50.000" if no sensor connected [°C].
- Env temp** Temperature from internal temperature sensor recording the environmental (room) temperature [°C].
- Peltier power** for regulation of block temperature, continuously recorded [% of maximum power].



7.4. Scaling

7.4.1. Graph \ Scaling F6

L[Ⓜ] Left mouse click on the X-, Y₁- or Y₂-axis opens the window **Graph layout Scaling**. **F6** provides flexibility to vary the

display of the plots and create Graph layouts. Viewing plots in differently scaled graphs, zooming the signal and time scales, and scrolling along the axes of the graph provide maximum information on the current experiment, but do not influence the format of stored data. Different ranges for the axes change the appearance of data dramatically.

Choose a Graph reference layout for a standard layout.

Select graph ▼ Pull down to select any of the defined graphs.

Minimum Defines the minimum axis position (suppression) for the display of data at a constant range.

Maximum Defines the maximum axis position for the display of data at a constant range.


Range Define the range of the Y_1 axis (left), Y_2 axis (right) and X axis (Time).

Select a predefined Graph layout from the options in the pull down menu. Edit the Graph layout name by a L[⌘] click on the name, and L[⌘] **Save** to save the entire graph layout.

Apply to Graph 1 The scaling defined for Graph 2 is applied to Graph 1.

7.4.2. Arrow keys

Use arrow keys to scroll (data on Y-axes), pan (time on X-axis), or zoom (expansion or compression: Ctrl+arrow key), independent of the **F6** window.



L[⌘] Select the active plot in a graph by a left click onto the label of the plot in the figure legend on the right. The active plot is highlighted. Scrolling and setting marks apply to the active plot.

↑ or ↓

Scroll up and down the Y axis, with a shift of 50% each time of the active plot.

Ctrl+↑

Magnify the signal (half the signal range is displayed).

Ctrl+↓

Cover a larger range on the screen for overview (twice the signal range is displayed).

→ or ←

Panning, shift 50% of the time axis to the right or left. During data acquisition, switch off automatic panning to pan backwards without changing the time range.

Ctrl+→


Magnify the time resolution on the screen (decreasing the time range). **Ctrl+→** expands the data, half the

original time range is displayed. The upper limit of the time range is fixed.

Ctrl+←

Cover a larger time range on the screen. **Ctrl+←** compresses the data, twice the original time range is displayed. The reference time point is fixed on the right during zooming in and out.

7.5. Graph \ Info



View and edit information, and load or save a graph layout (Section 8).

A quick selection of the tabs **Plots** **Ctrl+F6** and **Scaling** (**F6**) is possible.

7.6. Synchronous time axes

- Sets the time axes of all graphs at an identical range and offset, which is particularly useful while panning.

7.7. Autoscale time axis

Autoscale the entire experimental time scale.

7.8. Autoscale Y1 (Y2) axis

Autoscale the full data range.


7.9. Automatic pan

- Toggles automatic panning on/off, » [MiPNet19.18A](#), Status line.

7.10. Mouse control: Zoom **Ctrl+Z**

- Zoom in** Select **Mouse control: Zoom** in the menu or press **Ctrl+Z**. **L** click into the graph where you want to zoom in. Place the cursor at the upper lefthand corner of the field for zooming. Hold **Shift**, press the left mouse button **Shift+L** and slide the cursor to the lower righthand corner to define the field for zooming

in. The text for all axes now indicates the respective range (full scale).


-  **Zoom out** To return to the original scaling (as defined in the Graph layout), hold **Shift**, press the left mouse button and slide the cursor anywhere to the left **Shift+L**.

7.11. Mouse control: Mark **Ctrl+M**

- The Mark mode is active by default, or can be selected in the menu or by **Ctrl+M**. Specific sections of the experiment can be marked on each plot. Usually, marks are set on the plot for oxygen concentration for calibration ([MiPNet19.18D](#)), whereas marks on the plot for oxygen flux are set for exporting the median or average of flux to a table.

7.12. Full screen **On/Off**

- On** The selected graph may be shown alone on the full screen , or together with the other defined graphs . Full screen is particularly useful for a single channel overview and for **Copy to clipboard** **ALT+G B**.



7.13. Display numerical value

- On** The current numerical values are displayed in the graph for the active plots on the Y₁ axis and Y₂ axis (during data acquisition only).

7.14. Display Power-O2k

- On** The Power-O2k number, set in **Oxygraph-2k \ O2k configuration**, is shown in the active graph.


7.15. Graph options

- Show legends** **On/Off** Shows the quantities plotted on each Y axis.

Show gridlines gives options to show or hide vertical and horizontal gridlines.

- Show marks** **On/Off** optionally shows or hides marks in all graphs, without deleting the marks.

- Show events** **On/Off** optionally shows or hides events in all graphs, without deleting the events.



Font Change font size and style of axes labels and numbers.
Event line style can be modified.


- Show Oroboros icon in clipboard files** Yes!
- Show seconds on time axis** gives an option to exclude or include the display of seconds on the labels of the time axis.

7.16. Copy to clipboard

Select the active graph. In the Graph menu, **Copy to clipboard**, and select the WMF or BMP format (Graph\Copy to clipboard\WMF); **Alt+G B W**. Open the target file (DatLab-Excel template, Word, PowerPoint, Paint, etc.), and paste the image **Ctrl+V**.

9. Layout

Graph layouts are selected from the **Layout** menu for standardized display of graphs, plots and scaling of axes.



8.1. Info / Load / Save

Open this window (Section 7.5) to view and edit information, and load or save a graph layout.



Types

Four types of layout can be selected depending on the **O2k-configuration**:

- ▶ O2k-Core: Oxygen channel only.
- ▶ O2 & Amp: Oxygen and amperometric channel.
- ▶ O2 & pX: Oxygen and potentiometric channel.
- ▶ Other

8.2. Standard layouts – O2k-Core

In many Standard layouts, plots are shown for the left chamber in Graph 1 (top), and the right chamber in Graph 2 (below). During data acquisition, a 30 min time range is frequently used.

01 Calibration show Temp The layout used after switching on the O2k, for performing oxygen calibration. Oxygen concentration (blue lines, left Y-axis) and O2 slope negative (red lines, right Y-axis) are displayed in the top graph for the left chamber, and below for the right

chamber. Graph 3 (bottom) shows the block temperature on the left Y axis and the Peltier power on the right Y axis. The chambers reach thermal equilibrium when Peltier power stabilizes. Next observe equilibration of the oxygen signal with air in the gas phase above the stirred aqueous phase ('open' chamber), to perform an oxygen calibration ([MiPNet19.18D](#)).

02 Calibration - Background for recording O₂ sensor calibration and instrumental O₂ background test. 'O₂ slope neg.' is the negative slope of oxygen concentration, multiplied by 1000 to convert to units [pmol/ml], over time [s]. No correction is applied for instrumental O₂ background flux, $J^{\circ}_{O_2}$. 'O₂ slope neg.' is plotted on the right Y-axis with a scaling to display $\pm 10 \text{ pmol}\cdot\text{s}^{-1}\cdot\text{ml}^{-1}$ and zero in the middle of the Y₂ axis. Zero slope in the 'open' chamber at air calibration indicates stability of the oxygen signal. After closing the chamber, the slope deviates from zero as a function of the oxygen consumption of the polarographic oxygen sensor and of oxygen diffusion into or out of the chamber, which is the first point in the instrumental O₂ background test ([MiPNet19.18E](#)).

03 Background high O₂ for recording an instrumental O₂ background test at high oxygen from 150 to 450 μM .

04a Flux per volume displays background-corrected oxygen flux per volume, which is most relevant to evaluate experimental details, i.e. flux per volume is optimally in the range of 20 to 200-500 $\text{pmol}\cdot\text{s}^{-1}\cdot\text{ml}^{-1}$. This plot is also chosen when measurements on sample density are available only at a later stage ([MiPNet19.18E](#)). Total O₂ flux is corrected for instrumental O₂ background, $J^{\circ}_{O_2,V}$, to obtain sample oxygen flux per chamber volume, $J_{O_2,V}(\text{sample})$:

$$J_{O_2,V}(\text{sample}) = J_{O_2,V}(\text{total}) - J^{\circ}_{O_2,V}$$

04b Flux per volume overlay is similar to the layout '04a Flux per volume'. Graph 1 shows the background-corrected oxygen flux per volume [$\text{pmol}\cdot\text{s}^{-1}\cdot\text{ml}^{-1}$] superimposed from both chambers. Graph 2 shows oxygen concentration [μM].

05a Specific flux for plotting background-corrected oxygen flux per unit sample. The unit as a marker for the amount of sample is defined in the [F3](#) window. Example: Select Unit ▼ for the amount (mass) of sample added to

the chamber. The mass-concentration is automatically calculated (division by chamber volume, typically 2 ml). Then mass-specific flux [$\text{pmol}\cdot\text{s}^{-1}\cdot\text{mg}^{-1}$] is displayed as volume-specific flux [$\text{pmol}\cdot\text{s}^{-1}\cdot\text{ml}^{-1}$] divided by mass-concentration [mg/ml]:

$$J_{\text{O}_2,\text{mass}} = J_{\text{O}_2,\text{V}} / \text{mass-concentration}$$

05b Specific flux overlay is similar to '05a Specific flux'. Graph 1 shows the background-corrected oxygen flux per unit sample superimposed from both chambers. Graph 2 shows oxygen concentration [μM].

06a Specific flux high O2 is same as '05a Specific flux' with the scaling adapted to the high oxygen concentration used for permeabilized muscle fibers (150 to 450 μM ; range 300 μM).

06b Specific flux high O2 overlay as '06a Specific flux high O2', but oxygen flux (graph 1) and oxygen concentration (graph 2) superimposed from both chambers.

07a Flux control ratios shows the flux control ratio (*FCR*) and the oxygen concentration for the left chamber in graph 1 and for the right chamber in graph 2. First, the reference and baseline metabolic states are marked and the marks are named (Marks \ Names). Then the marks are selected in the menu Flux/Slope \ O2 slope.

07b Flux control ratios overlay as '07a Flux control ratios', but *FCR* (graph 1) and oxygen concentration (graph 2) superimposed from both chambers.

» Layouts for O2k-MultiSensor applications are explained in the specific sections of the O2k-Manual.

8.3. Lab layouts A Standard layout or any other layout can be modified and saved under Lab layouts, which is recommended for a team using project-specific layouts.

8.4. User: Name A Standard layout or any other layout can be modified and saved under a specific user name, which is recommended to distinguish individual layouts from standard layouts used by a team.

🔗 Further details:

http://www.bioblast.at/index.php/Select_plots_in_DatLab

10. Marks

Select

- Median
- Average
- Standard deviation
- Outlier index
- Range
- Maximum
- Minimum

Sort by

- Time
- Mark name

Set Marks to obtain the median, average, standard deviation, outlier index and range of the data within the mark, for calibration of the oxygen signal ([MiPNet19.18D](#)), flux analysis ([MiPNet19.18E](#)), or to delete marked data points. Marks define a section (period of time) in a selected plot (column) and contain the data of the selected plot within the defined section of time. Marks are shown by a horizontal bar in the active plot. Several marks can be set on any plot, but marks cannot overlap within a plot and are separated by one or more data points which are not marked.

Mark

CCP
Sample 1.1

- O2 concentration 1A
- O2 flow per cells(bc)1A

Shift+L

Shift+R

In the **Graph** menu select **Mouse control: Mark**, or press **Ctrl+M**.

Select the active plot in a graph by a **L** left click onto the label of the plot in the figure legend on the right. The active plot is highlighted.

Set the cursor on the starting position of the mark, hold **Shift** and press the left mouse button **Shift+L**, while moving the cursor along the X-axis.

The period of a mark may be reduced or the entire mark deleted by holding **Shift** and pressing the right mouse button **Shift+R**, while moving the cursor along the X-axis over the marked section.

Default mark names Marks are labelled with consecutive numbers by default, starting from 01, independent of the sequence of position of the marks. When deleting a mark **Shift+R**, then this mark number is missing on the mark list. Marks >99 are named as 00. Extending the default mark names, e.g. from '03' to '03-TD', maintains the numerical sequence for default names in subsequently generated marks.

Edit mark information

Start 01:13:43

Stop 01:16:03

N Points 70

Average 35.4613

Name 5P

Value 5.00000

Comment Pyruvate, 5 µl, 5 mM final

Delete points

Interpolate points

Recalc. slope

MitoPedia: Marks


Cancel OK

Mark information **L** click into the top or bottom bar of the mark:

Name: Edit the mark name.

Concentration: (*optional*) The value is used as the concentration or pX value in the Amp or pX calibration windows.

Volume: (*optional*) In SUIT protocols, the value is entered as the volume of a titration, for calculating the dilution of sample in DatLab-Excel templates.



Delete \ Interpolate points: Set marks for deleting disturbed flux during reoxygenations (Mark 03), or interpolating flux disturbed due to titrations (Marks 02 and 04; compare to Mark 1 without deletion or interpolation of points).
Restore points \ Recalc. flux: All data points can be restored (signal) or recalculated (flux).

9.1. Marks \ Statistics F2

Median	Unit	R	0Dig	1D	2Oct	3M.05	4c	5P	6G	10G
Value		0.000000	10.000000	10.000000	10.000000	2.000000	16.500000	5.000000	10.000000	10.000000
Start		00:10:53	00:30:41	00:37:46	00:42:16	00:47:08	01:10:27	01:13:43	01:18:07	01:20:10
Stop		00:13:58	00:32:40	00:39:48	00:43:42	00:48:25	01:11:32	01:16:03	01:20:10	
N Points		92	59	61	42	39	33	70	62	
6A: O2 concentration	µM	136.7183	102.4669	96.4319	92.1607	84.8638	118.7165	106.8785	93.8379	
6A: O2 slope neg.	pmol/(s*ml)	53.6480	15.1865	15.9293	24.3892	26.8240	51.4196	52.8227	52.6576	

1. Select channels for which plots should be displayed.
2. Select the O2k-chamber for which plots should be displayed in the marks statistics table.
3. Click on the source plot on which the marks are set.
4. **Traceability in DatLab-Excel templates:** Experimental details are copied to clipboard, where they can be edited further. Instrumental O₂ background correction and normalization of flux are calculated in the DatLab-Excel template, from O₂ slope neg. This allows for traceability of instrumental settings, normalization and baseline correction of flux in the spreadsheets. Instrumental background tests and analyses on the amount of sample can be completed after an



- experimental series, and updated corrections of flux can thus be finalized in the spreadsheet.
5. **Experimental details** are copied to clipboard. This option is de-selected for export of data as displayed in the graph windows, particularly for using DatLab-Excel templates without traceability.
 6. Select the statistical value calculated over the sections defined by marks. Default: **Median**.
 7. Select the sequence of marks, sorted as a **Time** sequence or **Mark name** in alphanumerical order.
 8. Values are displayed for selected plots. The source plot for marks is indicated by an "X".



L[⏏] click to select a single data cell of the table, copy by **Ctrl+C**, and paste into a Windows™ file by **Ctrl+V**.

Copy to clipboard: L[⏏] Copy to clipboard to copy the mark statistics table into an Excel or SigmaPlot file. It is important to carefully evaluate which set of data rows is relevant.


More details » [MiPNet19.18E O2 flux analysis](#).

Show system channels in statistics: Select or deselect system channels (barometric pressure, block temperature etc.) to be shown in Statistics.

9.2. Manage mark statistics setup

Rename or delete.

9.3. Specifications



Mark specifications (name, state, concentration and volume) are shown from the active plot. Define the **Template Name** and **Save template**.

Select a **Template**, L[⏏] click on **Preview from template**, and **Rename** all marks of the plot, including the state, concentration and volume.

9.4. Manage mark specifications templates

Rename or delete mark name templates in this window.

9.5. Copy marks from a selected plot to the active plot.



Updates » http://wiki.oroboros.at/index.php/MiPNet19.18C_DatLab-guide
 » http://wiki.oroboros.at/index.php/MiPNet19.18A_O2k-start

Next step – O2k-Core Manual D » [MiPNet19.18D O2k-calibration Supplement](#)

A. Some features of Windows™

- Print** The entire screen is 'printed' to clipboard.
- Alt+Print** The active window is 'printed' to clipboard.
- Ctrl+C** Copy to clipboard.
- Ctrl+V** The contents of the clipboard is pasted into the page of a Windows™ program (Word; Excel; PowerPoint).

B. Author contributions

- ¹Erich Gnaiger: ¹Oroboros Instruments; ²Universitätsklinik für Visceral-, Transplantations- und Thoraxchirurgie, D. Swarovski Forschungslabor, DSL, Medizinische Universität Innsbruck. EG contributed to the concept of all DatLab versions, wrote the chapter, edited the final version.
- ³Christina Plattner: Biozentrum Innsbruck - Sektion für Bioinformatik, SBI, Medizinische Universität Innsbruck. CP edited this chapter and complementary websites (help) on the ¹Oroboros -Bioblast homepage.
- ¹Ondrej Capek: ¹Oroboros Instruments. OC was mainly responsible for trouble shooting of test versions from DatLab 6 to DatLab 7 and contributed to writing this MiPNet publication.
- ⁴Lukas Gradl: software security networks – ssn, Innsbruck. LG is the programmer of all DatLab versions.

C. Acknowledgements



Quality Control in DatLab 7 is a contribution to the project MitoFit, funded by the Tyrolian Government within the programme K-Regio, Standortagentur Tirol, Austria.

<http://www.mitofit.org/index.php/O2k-MitoFit>

