



# Laboratory Protocol: Isolation of Rat Heart Mitochondria

Zuzana Sumbalova<sup>1</sup>, Mona Fontana-Ayoub<sup>2</sup>, Gerhard Krumschnabel<sup>2</sup>

<sup>1</sup> Pharmacobiochemical Laboratory of 3rd Department of Internal Medicine, Medical Faculty, Comenius University in Bratislava, Slovak Republic

<sup>2</sup>OROBOROS INSTRUMENTS Corp

high-resolution respirometry

Schöpfstr 18, A-6020 Innsbruck, Austria

Email: Gerhard.Krumschnabel@oroboros.at

[www.oroboros.at](http://www.oroboros.at)

**Preparation:** Switch on centrifuge and let it cool down to 4 °C. Keep all buffers, dissection gear, homogenization tools and centrifuge rotor at 4 °C (or on ice).

**Anesthesia:** Rats are anesthetised by intraperitoneal injection of Thiopental (0.1 g/kg) or by CO<sub>2</sub> narcosis.

## Isolation procedure:

1. kill rat, dissect out heart (take weight) and put it into ice-cold isolation medium A, wash to remove blood, discard all medium
2. cut the heart into small pieces (should become a mash), add drops of isolation buffer A while cutting
3. add isolation buffer B (10 ml/g tissue), transfer to pre-cooled glass/Teflon potter and homogenize with 8-12 strokes at medium speed (1000 rpm)
4. transfer the homogenate to a 50 ml beaker, bring volume up to ~ 15 ml/g tissue with isolation buffer A, place it on a magnetic stirrer, stir slowly for 20 min in an ice bath
5. re-homogenize the homogenate briefly in the potter, bring the volume up to ~ 20–30 ml with isolation buffer C
6. centrifuge: 1000 x g, 10 min, 4°C
7. transfer the supernatant to a new 50 ml falcon tube
8. centrifuge: 6200 x g, 10 min, 4°C
9. discard the supernatant, carefully re-suspend the mitochondrial pellet in a small volume of isolation buffer C and bring the volume to 20-30 ml (for 1 g tissue) with isolation buffer C
10. centrifuge: 6200 x g, 10 min, 4°C
11. discard supernatant and carefully re-suspend mitochondria with small volume of isolation buffer C (the volume of mitochondrial suspension for 1 g tissue ~ 1.2 ml)
12. store mitochondria in ice/water bath, use within 3-4 h
13. transfer subsamples (20 µl) into Eppendorf tubes and store at -20°C for further analysis (protein concentration, citrate synthase)

14. for respiration measurements add 2.5  $\mu$ l of mitochondrial suspension into a 2 ml chamber

**Isolation buffer A:**

Chemical	Final conc.	Required for 500 ml buffer
KCl	180 mM	6.71 g
EDTA	4 mM	0.745 g
BSA	1 g/l	0.5 g

Adjust pH to 7.4 with Tris, HCl

**Isolation buffer B:**

Isolation buffer A with 2.5 mg/ml Subtilisin. Add 25 mg Subtilisin to 10 ml of Buffer A.

**Isolation buffer C:**

Isolation buffer A without BSA

For 3 - 4 isolations 500 ml of media need to be prepared. Add BSA to 250 ml of medium and leave the residual 250 ml without BSA.